

## WORK BENCH: “WORKSTATION 1 AND 2”

- The workstation is a station with adjustable height to be used in the sterile packing area of sterilization facilities
- It can be static or mobile
- It is also a work bench used to control and verify the surgical instruments and the textile items and to pack the loads to be sterilized. It can be positioned between the washing and the sterilization area
- It can be completed with monitor, keyboard and mouse support in case you need to follow the traceability of the sterilization kits
- Height adjustability always provide a good ergonomic support for the operator’s best comfort
- Available in different sizes and versions
- The workstation can be customized with a wide range of equipment for the inspection and packing: roll dispenser with cutter, rolls holder, shelves, drawers, baskets, 220V sockets, data storage unit, automatic cut-out, etc

WORKSTATION



### WORKSTATION

| REF.          | WORK BENCH            |
|---------------|-----------------------|
| Workstation 1 | Static-without wheels |
| Workstation 2 | Mobile-with wheels    |

## WORKSTATION: “EC-DISPENSER” ROLL DISPENSER WITH CUTTER

- Roll dispenser with cutter for all types of rolls
- Suitable for either wall or table fastening
- Made of stainless steel

EC-DISPENSER



### EC-DISPENSER

| REF.              | SIZE (L x W x H)    | PACK  |
|-------------------|---------------------|-------|
| EC-DISPENSER 630  | 630 X 373 X 294 mm  | 1 pcs |
| EC-DISPENSER 910  | 910 X 373 X 294 mm  | 1 pcs |
| EC-DISPENSER 1260 | 1260 X 373 X 294 mm | 1 pcs |
| EC-DISPENSER 1510 | 1510 X 373 X 294 mm | 1 pcs |